For research use only

Version Number : 1.0-2110

# Foreasy Taq DNA Polymerase

HotStar Taq for PCR/qPCR

#### **Product Introduction**

Foreasy Taq DNA Polymerase is a new Taq enzyme expressed in Escherichia coli engineering bacteria by gene recombination technology. The enzyme itself has a certain hot-start activity and can be used for conventional PCR and qPCR; it has  $5'\rightarrow3'$  DNA polymerase activity and  $5'\rightarrow3'$  exonuclease activity, but no  $3'\rightarrow5'$  exonuclease activity.

# Component

Component	IM- 01011	IM-01012	IM-01013	
Foreasy Taq DNA Polymerase	5000 U (1 mL)	50 KU (10 mL)	500 KU (100 mL)	
(5 U/µL)				
2× Taq Reaction	25 mL ×5	250 mL ×5	×5 500 mL ×25	
Buffer				

#### Storage

-20  $\pm$  5 °C for 2 years or at -80 °C for long-term storage.

## Features

High specificity: The enzyme has a certain hot-start activity.

Fast Amplification: 10 sec/kb.

Highly adaptable template: can be used to efficiently amplify GC High value, various difficult-to-amplify DNA template.

Strong fidelity: ordinary Taq Enzyme 6 times.

Strong thermal stability: It can be placed at 37 °C for a week and maintains more than 90% activity

#### FOREGENE

# Application

Various PCR/qPCR systems and direct PCR systems PCR amplification of DNA fragments DNA labeling DNA sequencing PCR A-tailed

# **U** Definition

1U: The amount of enzyme required to incorporate 10 nmol of deoxynucleotides into acid-insoluble matter using activated salmon sperm DNA as template/primer for 30 minutes at 74°C.

## **Activity Assay Conditions**

1× Taq Reaction Buffer, 1.5 mM MgCl<sub>2</sub>; 0.2 mg/mL activated calf thymus DNA, 0.2 mM dNTPs $_{\circ}$ 

#### Storage

20 mM Tris-HCI (pH 8.0), 1 mM DTT, 0.1 mM EDTA, 50% glycerol, stabilizer。

## 2× Taq Reaction Buffer

Contains optimized ratios of Tris, KCI, MgCl2 and other ingredients.

## **Instance Reaction Protocol**

Template DNA	ΧμL
dNTPs (10 mM each)	1 µL
Primer-F	1 µL
Primer-R	1 µL
ddH2O	Το 50 μL
Total Volume	50 µL

#### **Reaction Condition**

Temperature	Time	Cycle	
37°C	5mins	1	
94°C	5mins	1	
94°C	10 Secs		
60°C	10 Secs	35	
72°C	20 sec/kb		
72°C	2mins	1	

**Note:** For 10  $\mu$ L and 20  $\mu$ L systems, add an equal volume of mineral oil if the thermal cycler does not have a heated lid.

PCR reaction conditions vary depending on the structural conditions of templates, primers, etc. In the specific operation ,it is necessary to design the optimal reaction conditions including the annealing temperature, extension time and so on , which is according to the different template type, the size of the target fragment, the base sequence of the amplified fragment, and the GC content and length of the primer.