

Mini Portable Centrifuge



FEATURES

- ◆ The machine material is high-strength engineering plastics.
- ◆ Brushless motor, ultra-low noise, stable operation.
- ◆ Digital Display. Closed loop control, real-time speed measurement, high speed accuracy.
- ◆ Adjustable acceleration time to adapt to different experiments.
- ◆ With electronic brake function, stop within 3 seconds, saving experiment time.
- ◆ Automatically open the door when stop.
- ◆ Time from 1s to 99min59s is adjustable.
- ◆ Speed 800rpm to the highest speed is adjustable.
- ◆ Adopt AC100-240V wide range power input.
- ◆ Three-in-one rotor (8*1.5/2ml+6*0.5ml+2*8*0.2ml), no need to replace the rotor.

TECHNICAL PARAMETER

Model	Mini-10K	Mini-7K	Mini-6k
Max Speed	10000r/min	7500r/min	6000r/min
Max RCF	6100xg	3450xg	2200xg
Rotor	Fixed angle rotor:8x1.5ml+6x0.5ml+2x8x0.2ml (Can fit flat tube and RIA tube)		
Noise	≤57dB(A)	≤52dB(A)	≤52dB(A)
Time Range	1s-99min59s		
Power Supply	AC 100-240V 50/60HZ		
Dimension	200x170x130mm		
Weight	1.2kg		